# **HD Freeform Basic Varifocals**

#### Product / Technical Information

- Introductory freeform lens with trusted traditional design
- Rx optimisation
- Balanced design between distance and near
- Non-compensated powers

# Ordering Requirements

Corridor length 16,18 Distance PD & Pupil Height


#### Patient Benefits

- Smooth transition between visual fields
- Quick and easy to adapt to design
- · Characteristics of traditional design with digital quality

# Material and Ranges

Clear		Photochromic Grey Brown		Polarised Grey Brown	
		Transitions Gen 8			
1.50	+7.00 to -8.00	1.50	+6.00 to -8.00	1.50	+6.00 to -8.00
1.60	+7.00 to -10.00	1.60	+6.00 to -9.00		
1.67	+8.00 to -10.00	1.67	+8.00 to -9.00		

Opposite cyls to 6.00 Adds from 1.00 to 3.00 in 0.25 steps



# HD Freeform Varifocals HD Standard, HD Defined, HD Personalised

#### Product / Technical Information

- Digital ray path technology reduces oblique astigmatism
- Steady lens design methodology
- Balanced design between distance and near

#### 3 Options:

HD Standard (prescribed power)

**HD Defined** (compensated power)

**HD Personalised** (individualised)

#### Patient Benefits

- Wide and balanced visual fields for all distances
- Quick and easy to adapt to design
- Smooth transition in all gaze directions
- Minimum lateral distortion
- Superior image stability

# Ordering Requirements

HD Standard, HD Defined: corridor length 14,16,18,20 distance PD & pupil height

HD Personalised: BVD, WRAP, PANTO Corridor length 14,16,18,20 distance PD & pupil height

# Material and Ranges

Clear		Transitions		Gen8	Xtractive
1.50	+9.00 to -10.00	1.50	+6.00 to -10.00	Brn Gry Grn Sphr Amthyst	Brn Gry Grn
1.53	+7.00 to -10.00	1.53	+7.00 to -9.00	Brn Gry	Brn Gry
1.59	+7.00 to -10.00	1.59	+7.00 to -11.00	Brn Gry	Brn Gy
1.60	+7.00 to -10.00	1.60	+7.00 to -10.00	Brn Gry Grn	Brn Gry
1.67	+10.00 to -12.00	1.67	+8.00 to -12.00	Brn Gry Grn	Brn Gry
1.74	+9.00 to -17.00	1.74	+10.00 to -15.00	Brn Gry	

Polari	Polarised						
1.50	+6.50 to -9.00	Brn Gry Grn, Drivewear, Xtractive polarised					
1.53	+7.00 to -6.00	Brn Gry					
1.59	+7.00 to -8.00	Brn Gry Grn, Drivewear, Xtractive polarised, Graduated Brn Gry					
1.60	+8.00 to -10.00	Brn Grn Gry, Xtractive polarised					
1.67	+8.00 to -12.00	Brn Grn Gry, Xtractive polarised					

Opposite cyls to 6.00 Adds from 0.50 to 4.50 in 0.25 steps



# **Coatings Available**

- Kleencoat plus (2 year scratch warranty)
- Kleencoat
- Kleencoat Blue
- Kleencoat Sun
- Kleencoat Drive

## **Mirror Coat**

Silver, Gold, Blue, Fire, Pink, Green (1.50, 1.60, 1.67 index only)

## **Blue Protect System**



# HD Freeform Camber Varifocals HD Camber Defined, HD Camber Personalised

#### Product / Technical Information

- Camber variable front curve increasing from top to bottom of the lens
- Individual frame fit and wearers preferences
- Digital ray path technology reduces oblique astigmatism
- Steady lens design methodology
- Balanced soft lens design

#### 2 Options:

HD Camber Defined (compensated)
HD Camber Personalised (individualised)

#### Patient Benefits

- Significantly improved reading and arm's length vision fields
- · Smooth transition in all gaze directions
- Minimum lateral distortion
- Superior image stability

## Ordering Requirements

HD Camber Defined: Corridor length 14,16,18,20 Distance PD & pupil height

HD Camber Personalised: BVD, WRAP, PANTO Corridor length 14,16,18,20 Distance PD & pupil height

# Material and Ranges

Clear	Clear		Transitions Gen8 Grey Brown		Polarised Grey Brown	
1.50	+6.00 to -11.00	1.50	+6.00 to -11.00	1.50	+6.00 to -11.00	
1.53	+6.00 to -10.00	1.53	+6.00 to -10.00			
1.59	+6.00 to -10.00	1.59	+6.00 to -10.00			
1.60	+7.00 to -11.00	1.60	+7.00 to -11.00	1.60	+7.00 to -11.00	
1.67	+7.00 to -11.00	1.67	+7.00 to -11.00	1.67	+7.00 to -11.00	
1.74	+8.00 to -12.00	1.74	+8.00 to -12.00			

Opposite cyls to 4.00

Adds from 0.50 to 4.50 in 0.25 steps



# **Coatings Available**

- Kleencoat plus (2 year scratch warranty)
- Kleencoat
- Kleencoat Blue
- Kleencoat Sun
- Kleencoat Drive

#### **Mirror Coat**

Silver, Gold, Blue, Fire, Pink, Green (1.50, 1.60, 1.67 index only)

## **Blue Protect System**



# **HD Camber Life Varifocal**

#### Product / Technical Information

- Patient lifestyle information used to customise design
- Camber variable front curve increasing from top to bottom of the lens
- Individual frame fit and wearers preferences
- Digital ray path technology reduces oblique astigmatism
- Steady lens design methodology
- Extra soft lens design

#### Patient Benefits

- Custom designed to meet patient's lifestyle requirements
- Significantly improved reading and arm's length vision fields
- Minimum lateral distortion
- Superior image stability, less swim / sway

# Ordering Requirements

- State prolonged use for distance, intermediate, reading, first time or experienced wearer
- BVD, Panto, Wrap
- Corridor length 14,16,18,20mm
- Distance PD & pupil height

# Material and Ranges

Clear		Transitions Grey Brown		Polarised Grey Brown	
1.50	+6.00 to -11.00	1.50	+6.00 to -11.00	1.50	+6.00 to -11.00
1.53	+6.00 to -10.00	1.53	+6.00 to -10.00		
1.59	+6.00 to -10.00	1.59	+6.00 to -10.00		
1.60	+7.00 to -11.00	1.60	+7.00 to -11.00	1.60	+7.00 to -11.00
1.67	+7.00 to -11.00	1.67	+7.00 to -11.00	1.67	+7.00 to -11.00
1.74	+8.00 to -12.00	1.74	+8.00 to -12.00		

Opposite cyls to 4.00 Adds from 0.50 to 4.50 in 0.25 steps



# **Coatings Available**

- Kleencoat plus (2 year scratch warranty)
- Kleencoat
- Kleencoat Blue
- Kleencoat Sun
- Kleencoat Drive

#### **Mirror Coat**

Silver, Gold, Blue, Fire, Pink, Green (1.50, 1.60, 1.67 index only)

# **Blue Protect System**



# **HD Sun8**

# **Single Vision & Varifocal**

#### Product / Technical Information

- Front curve up to base 8
- Oblique astigmatism correction for expanded peripheral vision, ideal for large frames
- Optimised freeform centration capable of producing blank diameters equivalent to 110mm, ideal for large frames
- Digital ray path technology

#### Patient Benefits

- Wide fields of view
- Expertly produced lenses that maintain the original look & design of the original plano sun spec

# Ordering Requirements

PDs and pupil heights


# Material and Ranges

Clear		Transitions		Gen8	Xtractive
1.50	+6.00 to -6.00	1.50	+6.00 to -6.00	Brn Gry Grn Saph Amthyst	Brn Gry Grn
1.59	+6.00 to -6.00	1.59	+6.00 to -6.00	Brn Gry	Brn Gry
1.60	+6.00 to -6.00	1.60	+6.00 to -6.00	Brn Gry Grn	Brn Gry
1.67	+6.00 to -6.00	1.67	+6.00 to -6.00	Brn Gry Grn	Brn Gry

Polarised						
1.50	+6.00 to -6.00	Brn Gry Grn, Drivewear, Xtractive Polarised				
1.59	+6.00 to -6.00	Brn Gry Grn, Drivewear, Xtractive Polarised, Graduated Brn Gry				
1.60	+6.00 to -6.00	Brn Grn Gry, Xtractive Polarised				
1.67	+6.00 to -6.00	Brn Grn Gry, Xtractive Polarised				

Opposite cyls to 4.0

Adds from 1.00 to 3.50 in 0.25 steps





# **Coatings Available**

- Kleencoat Sun
- Kleencoat Drive

## **Mirror Coat**

Silver, Gold, Blue, Fire, Pink Gold, Green, Orange

# **HD Drive**

# Single Vision & Varifocal

# Product / Technical Information

- Optimised distance and intermediate zones
- Night zone, helps compensate the refractive change that occurs between day and night
- Blue Protect System highly recommended for night time glare Digital ray path technology

#### Patient Benefits

- Wider distance and intermediate visual fields
- Better vision of the dashboard and external mirrors
- Reduced fatigue when night driving

# Ordering Requirements

PDs and pupil heights

# Material and Ranges

Clear		Trans	itions	Gen8	Xtractive
1.50	+9.00 to -10.00	1.50	+6.00 to -10.00	Brn Gry Grn Saph Amthyst	Brn Gry Grn
1.53	+7.00 to -10.00	1.53	+7.00 to -9.00	Brn Gry	Brn Gry
1.59	+7.00 to -10.00	1.59	+7.00 to -11.00	Brn Gry	Brn Gry
1.60	+7.00 to -10.00	1.60	+7.00 to -10.00	Brn Gry Grn	Brn Gry
1.67	+10.00 to -12.00	1.67	+8.00 to -12.00	Brn Gry Grn	Brn Gry
1.74	+9.00 to -17.00	1.74	+10.00 to -15.00	Brn Gry	

Polari	Polarised							
1.50	+6.50 to -9.00	Brn Gry Grn	Drivewear					
1.53	+7.00 to -6.00	Brn Gry						
1.59	+7.00 to -8.00	Brn Gry Grn,	Drivewear					
1.60	+8.00 to -10.00	Brn Grn Gry						
1.67	+8.00 to -12.00	Brn Grn Gry						

Opposite cyls to 6.0

Adds from 0.50 to 4.50 in 0.25 steps



# **Coatings Available**

- Kleencoat plus (2 year scratch warranty)
- Kleencoat Blue
- Kleencoat Sun
- Kleencoat Drive

#### Mirror Coat

Silver, Gold, Blue, Fire, Pink, Green (1.50, 1.60, 1.67 index only)

# **Blue Protect System**

